

For Research Use Only

RTCB Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02340-2

Capture Antibody Information

Catalog Number:
86258-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250866B12
Reactivity:
human
Immunogen Catalog Number:
Ag13836

Conjugate:
Unconjugated
Full name:
chromosome 22 open reading frame 28
Gene ID:
51493

Detection Antibody Information

Catalog Number:
86258-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250866H8
Reactivity:
human
GenBank:
BC000151
Immunogen Catalog Number:
Ag13836

Conjugate:
Unconjugated
Full name:
chromosome 22 open reading frame 28
Gene ID:
51493

Applications

Tested Applications:
Sandwich ELISA

Range:
0.391-25 ng/ml (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02340-2 targets RTCB in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: C22orf28 Recombinant antibody, PBS Only (Capture) 86258-3-PBS (250866B12). 100 µg. Concentration 1 mg/mL.

Detection antibody: C22orf28 Recombinant antibody, PBS Only (Detector) 86258-1-PBS (250866H8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02340-2, Human RTCB Recombinant Matched Antibody Pair - PBS only. 86258-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag13836. 86258-1-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean RTCB concentration was determined to be 78.30 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load and 77.70 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.