

# Alpha Actinin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02360-1**

## Capture Antibody Information

Catalog Number:  
86340-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
250833F4

Reactivity:  
human

Immunogen Catalog Number:  
Ag1859

Conjugate:  
Unconjugated

Full name:  
actinin, alpha 1

Gene ID:  
87

## Detection Antibody Information

Catalog Number:  
86340-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
250833B7

Reactivity:  
human

GenBank:  
BC015766

Immunogen Catalog Number:  
Ag1859

Conjugate:  
Unconjugated

Full name:  
actinin, alpha 1

Gene ID:  
87

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02360-1 targets Alpha Actinin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: a-Actinin Recombinant antibody, PBS Only (Capture) 86340-2-PBS (250833F4). 100 µg. Concentration 1 mg/mL.

Detection antibody: a-Actinin Recombinant antibody, PBS Only (Detector) 86340-1-PBS (250833B7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Alpha Actinin concentration was determined to be 254.6 ng/mL in NIH/3T3 cell extract based on a 1.4 mg/mL extract load.



Cytometric bead array standard curve of MP02360-1, Alpha Actinin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86340-2-PBS. Detection antibody: 86340-1-PBS. Standard: Ag1859. Range: 0.781-100 ng/mL.