

CSNK1D Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02443-1**

Capture Antibody Information

Catalog Number: 86423-3-PBS	Clone ID: 251238B11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: casein kinase 1, delta
Isotype: IgG	Immunogen Catalog Number: Ag5856	Gene ID: 1453
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 86423-1-PBS	Clone ID: 251238G5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: casein kinase 1, delta
Isotype: IgG	GenBank: BC003558	Gene ID: 1453
Purification Method: Protein A purification	Immunogen Catalog Number: Ag5856	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP02443-1 targets CSNK1D in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CSNK1D Recombinant monoclonal antibody, PBS Only (Capture) 86423-3-PBS (251238B11). 100 µg. Concentration 1 mg/mL.

Detection antibody: Casein Kinase 1 Delta Recombinant monoclonal antibody, PBS Only (Detector) 86423-1-PBS (251238G5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02443-1, CSNK1D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86423-3-PBS. Detection antibody: 86423-1-PBS. Standard: Ag5856. Range: 0.625-80 ng/mL



The mean CSNK1D concentration was determined to be 12.5 ng/mL in HEK-293 cell extract based on a 1.1 mg/mL extract load.