

For Research Use Only



# Zyxin Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP02506-1

## Capture Antibody Information

Catalog Number:  
86364-5-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
250924H4

Reactivity:  
human

Immunogen Catalog Number:  
Ag0419

Conjugate:  
Unconjugated

Full name:  
zyxin

Gene ID:  
7791

## Detection Antibody Information

Catalog Number:  
86364-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
250924C7

Reactivity:  
human

GenBank:  
BC008743

Immunogen Catalog Number:  
Ag0419

Conjugate:  
Unconjugated

Full name:  
zyxin

Gene ID:  
7791

## Applications

Tested Applications:  
Sandwich ELISA

Range:  
0.625-40 ng/mL (Sandwich ELISA)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02506-1 targets Zyxin in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: ZYX Recombinant monoclonal antibody, PBS Only (Capture) 86364-5-PBS (250924H4). 100 µg. Concentration 1 mg/mL.

Detection antibody: ZYX Recombinant monoclonal antibody, PBS Only (Detector) 86364-4-PBS (250924C7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02506-1, Human Zyxin Recombinant Matched Antibody Pair - PBS only. 86364-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0419. 86364-4-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL



Plasma of six individual healthy human donors was measured. The human Zyxin concentration of detected samples was determined to be 69.10 ng/mL with a range of 25.02 - 134.21 ng/mL