

Mouse Cd2ap Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02649-1**

Capture Antibody Information

Catalog Number:
86727-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251428H6

Reactivity:
mouse

Immunogen Catalog Number:
Ag0622

Conjugate:
Unconjugated

Full name:
CD2-associated protein

Gene ID:
12488

Detection Antibody Information

Catalog Number:
86727-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251428G6

Reactivity:
mouse

GenBank:
BC157901

Immunogen Catalog Number:
Ag0622

Conjugate:
Unconjugated

Full name:
CD2-associated protein

Gene ID:
12488

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02649-1 targets Cd2ap in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Cd2ap Recombinant antibody, PBS Only (Capture) 86727-3-PBS (251428H6). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Cd2ap Recombinant antibody, PBS Only (Detector) 86727-2-PBS (251428G6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Cd2ap concentration was determined to be 948.7 ng/mL in Mouse kidney tissue extract based on a 8.4 mg/mL extract load.



Cytometric bead array standard curve of MP02649-1, MOUSE Cd2ap Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86727-3-PBS. Detection antibody: 86727-2-PBS. Standard: Ag0622. Range: 0.781-100 ng/mL.