

Entactin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02707-1**

Capture Antibody Information

Catalog Number:

86702-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

251592D3

Reactivity:

human

Immunogen Catalog Number:

Eg3912

Conjugate:

Unconjugated

Full name:

nidogen 1

Gene ID:

4811

Detection Antibody Information

Catalog Number:

86702-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

251592C11

Reactivity:

human

GenBank:

NM_002508.3

Immunogen Catalog Number:

Eg3912

Conjugate:

Unconjugated

Full name:

nidogen 1

Gene ID:

4811

Applications

Tested Applications:

Sandwich ELISA

Range:

1.56-50 ng/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02707-1 targets Entactin in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: NID1 Recombinant monoclonal antibody, PBS Only (Capture) 86702-3-PBS (251592D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: NID1 Recombinant monoclonal antibody, PBS Only (Detector) 86702-2-PBS (251592C11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02707-1, Human Entactin Recombinant Matched Antibody Pair - PBS only. 86702-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3912. 86702-2-PBS was HRP conjugated as the detection antibody. Range: 1.56-50 ng/mL.



Plasma of six individual healthy human donors was measured. The Entactin concentration of detected samples was determined to be 32.17 ng/mL with a range of 25.81-44.92 ng/mL.