

For Research Use Only

# UNG Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02750-1

## Capture Antibody Information

Catalog Number:  
86880-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
251858H8

Reactivity:  
human

Immunogen Catalog Number:  
Ag3072

Conjugate:  
Unconjugated

Full name:  
uracil-DNA glycosylase

Gene ID:  
7374

## Detection Antibody Information

Catalog Number:  
86880-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
251858E9

Reactivity:  
human

GenBank:  
BC015205

Immunogen Catalog Number:  
Ag3072

Conjugate:  
Unconjugated

Full name:  
uracil-DNA glycosylase

Gene ID:  
7374

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.562-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02750-1 targets UNG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: UNG Recombinant monoclonal antibody, PBS Only (Capture) 86880-2-PBS (251858H8). 100 µg. Concentration 1 mg/mL.

Detection antibody: UNG Recombinant monoclonal antibody, PBS Only (Detector) 86880-1-PBS (251858E9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MPO2750-1, UNG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86880-2-PBS. Detection antibody: 86880-1-PBS. Standard: Ag3072. Range: 1.562-200 ng/mL.



The mean UNG concentration was determined to be 3284.0 ng/mL in HT-29 cell extract based on a 1.5 mg/mL extract load, 1758.9 ng/mL in HuH-7 cell extract based on a 1.2 mg/mL extract load, 2609.9 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.