

For Research Use Only

# SLC27A3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02758-1

## Capture Antibody Information

**Catalog Number:**  
86047-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250628G3  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag4015

**Conjugate:**  
Unconjugated  
**Full name:**  
solute carrier family 27 (fatty acid transporter), member 3  
**Gene ID:**  
11000

## Detection Antibody Information

**Catalog Number:**  
86047-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250628F10  
**Reactivity:**  
human  
**GenBank:**  
BC029792  
**Immunogen Catalog Number:**  
Ag4015

**Conjugate:**  
Unconjugated  
**Full name:**  
solute carrier family 27 (fatty acid transporter), member 3  
**Gene ID:**  
11000

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.312-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02758-1 targets SLC27A3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC27A3 Recombinant monoclonal antibody, PBS Only (Capture) 86047-3-PBS (250628G3). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC27A3 Recombinant monoclonal antibody, PBS Only (Detector) 86047-2-PBS (250628F10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.



## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02758-1, SLC27A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86047-3-PBS. Detection antibody: 86047-2-PBS. Standard: Ag4015. Range: 0.312-40 ng/mL.

The mean SLC27A3 concentration was determined to be 1.1 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.