

For Research Use Only

EIF2C1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02777-1

Capture Antibody Information

Catalog Number:
80568-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251638D7
Reactivity:
human

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 2C, 1
Gene ID:
26523

Detection Antibody Information

Catalog Number:
80568-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251638C8
Reactivity:
human
GenBank:
NM_012199

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 2C, 1
Gene ID:
26523

Applications

Tested Applications:
Sandwich ELISA

Range:
78.1-5000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02777-1 targets EIF2C1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: EIF2C1-Specific Recombinant monoclonal antibody, PBS Only (Capture) 80568-3-PBS (251638D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: EIF2C1-Specific Recombinant monoclonal antibody, PBS Only (Detector) 80568-2-PBS (251638C8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02777-1, Human EIF2C1 Recombinant Matched Antibody Pair - PBS only. 80568-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00435. 80568-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



The mean EIF2C1 concentration was determined to be 111.5 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.