

Vitronectin/S-protein Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02822-1**

Capture Antibody Information

Catalog Number:

87005-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252166C10

Reactivity:

human

Immunogen Catalog Number:

Eg4682

Conjugate:

Unconjugated

Full name:

vitronectin

Gene ID:

7448

Detection Antibody Information

Catalog Number:

87005-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252166H1

Reactivity:

human

GenBank:

NM_000638.4

Immunogen Catalog Number:

Eg4682

Conjugate:

Unconjugated

Full name:

vitronectin

Gene ID:

7448

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02822-1 targets Vitronectin/S-protein in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Vitronectin Recombinant monoclonal antibody, PBS Only (Capture) 87005-2-PBS (252166C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: Vitronectin Recombinant monoclonal antibody, PBS Only (Detector) 87005-1-PBS (252166H1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02822-1, Vitronectin/S-protein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87005-2-PBS. Detection antibody: 87005-1-PBS. Standard: Eg4682. Range: 0.625-80 ng/mL



The mean Vitronectin/S-protein concentration was determined to be 430.3 ug/mL in human plasma.