

For Research Use Only

NUSAP1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02838-1

Capture Antibody Information

Catalog Number:
87034-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252131D5
Reactivity:
human
Immunogen Catalog Number:
Ag2654

Conjugate:
Unconjugated
Full name:
nucleolar and spindle associated protein 1
Gene ID:
51203

Detection Antibody Information

Catalog Number:
87034-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252131D10
Reactivity:
human
GenBank:
BC001308
Immunogen Catalog Number:
Ag2654

Conjugate:
Unconjugated
Full name:
nucleolar and spindle associated protein 1
Gene ID:
51203

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02838-1 targets NUSAP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NUSAP1 Recombinant monoclonal antibody, PBS Only (Capture) 87034-2-PBS (252131D5). 100 µg. Concentration 1 mg/mL.

Detection antibody: NUSAP1 Recombinant monoclonal antibody, PBS Only (Detector) 87034-1-PBS (252131D10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.



Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2838-1, NUSAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87034-2-PBS. Detection antibody: 87034-1-PBS. Standard: Ag2654. Range: 0.625-80 ng/mL.

The mean NUSAP1 concentration was determined to be 99.0 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 84.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 38.4 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.