

CNPase Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02922-1**

Capture Antibody Information

Catalog Number:
86967-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252105A6
Reactivity:
human
Immunogen Catalog Number:
Ag4227

Conjugate:
Unconjugated
Full name:
2',3'-cyclic nucleotide 3'
phosphodiesterase
Gene ID:
1267

Detection Antibody Information

Catalog Number:
86967-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252105D4
Reactivity:
human
GenBank:
BC028040
Immunogen Catalog Number:
Ag4227

Conjugate:
Unconjugated
Full name:
2',3'-cyclic nucleotide 3'
phosphodiesterase
Gene ID:
1267

Applications

Tested Applications:
Sandwich ELISA

Range:
7.81-500 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02922-1 targets CNPase in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: CNP Recombinant monoclonal antibody, PBS Only (Capture) 86967-6-PBS (252105A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: CNP Recombinant monoclonal antibody, PBS Only (Detector) 86967-3-PBS (252105D4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02922-1, Human CNPase Recombinant Matched Antibody Pair - PBS only. 86967-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag4227. 86967-3-PBS was HRP conjugated as the detection antibody. Range: 7.81-500 ng/mL



The mean CNPase concentration was determined to be 3,335.70 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load and 2,058.70 ng/mL in HCT 116 cell extract based on a 1.30 mg/mL extract load.