

For Research Use Only

RCOR1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02924-1

Capture Antibody Information

Catalog Number:
86954-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252115G5

Reactivity:
human

Immunogen Catalog Number:
Ag26770

Conjugate:
Unconjugated

Full name:
REST corepressor 1

Gene ID:
23186

Detection Antibody Information

Catalog Number:
86954-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252115F6

Reactivity:
human

GenBank:
NM_015156

Immunogen Catalog Number:
Ag26770

Conjugate:
Unconjugated

Full name:
REST corepressor 1

Gene ID:
23186

Applications

Tested Applications:
Sandwich ELISA

Range:
0.781-50 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02924-1 targets RCOR1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: RCOR1 Recombinant monoclonal antibody, PBS Only (Capture) 86954-3-PBS (252115G5). 100 µg. Concentration 1 mg/mL.

Detection antibody: RCOR1 Recombinant monoclonal antibody, PBS Only (Detector) 86954-2-PBS (252115F6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02924-1, Human RCOR1 Recombinant Matched Antibody Pair - PBS only. 86954-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag26770. 86954-2-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL



The mean RCOR1 concentration was determined to be 903.00 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load, 279.00 ng/mL in K-562 cell extract based on a 1.20 mg/mL extract load, 566.90 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load and 267.70 ng/mL in HuH-7 cell extract based on a 1.20 mg/mL extract load.