

For Research Use Only

IRF4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02975-1

Capture Antibody Information

Catalog Number:

87342-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252430B2

Reactivity:

human

Immunogen Catalog Number:

Ag1762

Conjugate:

Unconjugated

Full name:

IRF 4

Gene ID:

3662

Detection Antibody Information

Catalog Number:

87342-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252430D12

Reactivity:

human

GenBank:

BC015752

Immunogen Catalog Number:

Ag1762

Conjugate:

Unconjugated

Full name:

IRF 4

Gene ID:

3662

Applications

Tested Applications:

Cytometric bead array

Range:

0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02975-1 targets IRF4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IRF4 Recombinant monoclonal antibody, PBS Only (Capture) 87342-1-PBS (252430B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: IRF4 Recombinant monoclonal antibody, PBS Only (Detector) 87342-2-PBS (252430D12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02975-1, IRF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87342-1-PBS. Detection antibody: 87342-2-PBS. Standard: Ag1762. Range: 0.625-80 ng/mL

The mean IRF4 concentration was determined to be 58.7 ng/mL in Ramos cell extract based on a 2.3 mg/mL extract load, THP-1 cell extract were served as a negative control.