

For Research Use Only



POFUT1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP03009-1

Capture Antibody Information

Catalog Number:
87434-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252682C5

Reactivity:
human

Immunogen Catalog Number:
Eg6687

Conjugate:
Unconjugated

Full name:
protein O-fucosyltransferase 1

Gene ID:
23509

Detection Antibody Information

Catalog Number:
87434-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252682G5

Reactivity:
human

GenBank:
NM_015352.2

Immunogen Catalog Number:
Eg6687

Conjugate:
Unconjugated

Full name:
protein O-fucosyltransferase 1

Gene ID:
23509

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03009-1 targets POFUT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: POFUT1/FUT12 Recombinant monoclonal antibody, PBS Only (Capture) 87434-2-PBS (252682C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: POFUT1 Recombinant monoclonal antibody, PBS Only (Detector) 87434-1-PBS (252682G5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03009-1, POFUT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87434-2-PBS. Detection antibody: 87434-1-PBS. Standard: Eg6687. Range: 0.625-80 ng/mL.

The mean POFUT1 concentration was determined to be 89.9 ng/mL in CaCO-2 cell extract based on a 1.3 mg/mL extract load, 9.2 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.