

For Research Use Only

# HGS Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03019-1

## Capture Antibody Information

**Catalog Number:**  
87420-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
252813A5  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag0589

**Conjugate:**  
Unconjugated  
**Full name:**  
hepatocyte growth factor-regulated tyrosine kinase substrate  
**Gene ID:**  
9146

## Detection Antibody Information

**Catalog Number:**  
87420-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
252813F1  
**Reactivity:**  
human  
**GenBank:**  
BC003565  
**Immunogen Catalog Number:**  
Ag0589

**Conjugate:**  
Unconjugated  
**Full name:**  
hepatocyte growth factor-regulated tyrosine kinase substrate  
**Gene ID:**  
9146

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP03019-1 targets HGS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HGS Recombinant monoclonal antibody, PBS Only (Capture) 87420-3-PBS (252813A5). 100 µg. Concentration 1 mg/mL.

Detection antibody: HGS Recombinant monoclonal antibody, PBS Only (Detector) 87420-2-PBS (252813F1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP03019-1, HGS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87420-3-PBS. Detection antibody: 87420-2-PBS. Standard: Ag0589. Range: 0.781-100 ng/mL.



The mean HGS concentration was determined to be 346.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 327.2 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load, mouse skeletal muscle tissue extract were served as a negative control.