

NCOA6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03030-2**

Capture Antibody Information

Catalog Number:
87454-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252826E2

Reactivity:
human

Immunogen Catalog Number:
Ag18737

Conjugate:
Unconjugated

Full name:
nuclear receptor coactivator 6

Gene ID:
23054

Detection Antibody Information

Catalog Number:
87454-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252826G11

Reactivity:
human

GenBank:
BC136272

Immunogen Catalog Number:
Ag18737

Conjugate:
Unconjugated

Full name:
nuclear receptor coactivator 6

Gene ID:
23054

Applications

Tested Applications:
Sandwich ELISA

Range:
19.5-1250 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03030-2 targets NCOA6 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: NCOA6 Recombinant monoclonal antibody, PBS Only (Capture) 87454-3-PBS (252826E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: NCOA6 Recombinant monoclonal antibody, PBS Only (Detector) 87454-4-PBS (252826G11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP03030-2, Human NCOA6 Recombinant Matched Antibody Pair - PBS only. 87454-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag18737. 87454-4-PBS was HRP conjugated as the detection antibody. Range: 19.5-1250 pg/mL



The mean NCOA6 concentration was determined to be 8.9 ng/mL in BxPC-3 cell extract based on a 1.2 mg/mL extract load, 12.0 ng/mL in HT-29 cell extract based on a 1.1 mg/mL extract load and 6.4 ng/mL in Romas cell extract based on a 1.2 mg/mL extract load.