

For Research Use Only

# Mouse CDCP1/CD318 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03057-1

## Capture Antibody Information

Catalog Number:  
87514-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
251667A2

Reactivity:  
mouse

Immunogen Catalog Number:  
Eg6199

Conjugate:  
Unconjugated

Full name:  
CUB domain containing protein 1

Gene ID:  
109332

## Detection Antibody Information

Catalog Number:  
87514-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
251667G9

Reactivity:  
mouse

GenBank:  
NM\_133974.3

Immunogen Catalog Number:  
Eg6199

Conjugate:  
Unconjugated

Full name:  
CUB domain containing protein 1

Gene ID:  
109332

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.312-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP03057-1 targets CDCP1/CD318 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD318/CDCP1 Recombinant monoclonal antibody, PBS Only (Capture) 87514-3-PBS (251667A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD318/CDCP1 Recombinant monoclonal antibody, PBS Only (Detector) 87514-2-PBS (251667G9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP03057-1, MOUSE CDCP1/CD318 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87514-3-PBS. Detection antibody: 87514-2-PBS. Standard: Eg6199. Range: 0.312-40 ng/mL



The mean CDCP1/CD318 concentration was determined to be 13.9 ng/mL in mouse skin tissue extract based on a 3.9 mg/mL extract load, mouse skeletal muscle tissue extract were served as a negative control.