

For Research Use Only

BRWD1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03069-1

Capture Antibody Information

Catalog Number:
87591-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252838E5

Reactivity:
human

Immunogen Catalog Number:
Ag30902

Conjugate:
Unconjugated

Full name:
bromodomain and WD repeat domain containing 1

Gene ID:
54014

Detection Antibody Information

Catalog Number:
87591-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
252838H1

Reactivity:
human

GenBank:
BC064602

Immunogen Catalog Number:
Ag30902

Conjugate:
Unconjugated

Full name:
bromodomain and WD repeat domain containing 1

Gene ID:
54014

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03069-1 targets BRWD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BRWD1 Recombinant monoclonal antibody, PBS Only (Capture) 87591-2-PBS (252838E5). 100 µg. Concentration 1 mg/mL.

Detection antibody: BRWD1 Recombinant monoclonal antibody, PBS Only (Detector) 87591-1-PBS (252838H1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03069-1, BRWD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87591-2-PBS. Detection antibody: 87591-1-PBS. Standard: Ag30902. Range: 0.156-20 ng/mL

The mean BRWD1 concentration was determined to be 4.6 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 8.6 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 2.5 ng/mL in U2OS cell extract based on a 1.2 mg/mL extract load.