

For Research Use Only

# RIC8B Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03094-1

## Capture Antibody Information

Catalog Number:

87524-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252908B11

Reactivity:

human

Immunogen Catalog Number:

Ag11888

Conjugate:

Unconjugated

Full name:

resistance to inhibitors of cholinesterase 8 homolog B (C. elegans)

Gene ID:

55188

## Detection Antibody Information

Catalog Number:

87524-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252908A7

Reactivity:

human, mouse, rat, pig

GenBank:

BC034220

Immunogen Catalog Number:

Ag11888

Conjugate:

Unconjugated

Full name:

resistance to inhibitors of cholinesterase 8 homolog B (C. elegans)

Gene ID:

55188

## Applications

Tested Applications:

Cytometric bead array

Range:

0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP03094-1 targets RIC8B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RIC8B Recombinant monoclonal antibody, PBS Only (Capture) 87524-2-PBS (252908B11). 100 µg. Concentration 1 mg/mL.

Detection antibody: RIC8B Recombinant monoclonal antibody, PBS Only (Detector) 87524-1-PBS (252908A7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:



T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP03094-1, RIC8B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87524-2-PBS. Detection antibody: 87524-1-PBS. Standard: Ag11888. Range: 0.156-20 ng/mL

The mean RIC8B concentration was determined to be 1.4 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 3.9 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.