

For Research Use Only

IL-17C Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03144-1

Capture Antibody Information

Catalog Number:
87706-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251249G4

Reactivity:
human

Immunogen Catalog Number:
Eg2788

Conjugate:
Unconjugated

Full name:
interleukin 17C

Gene ID:
27189

Detection Antibody Information

Catalog Number:
87706-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251249F9

Reactivity:
human

GenBank:
NM_013278.3

Immunogen Catalog Number:
Eg2788

Conjugate:
Unconjugated

Full name:
interleukin 17C

Gene ID:
27189

Applications

Tested Applications:
Sandwich ELISA

Range:
15.6-1000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03144-1 targets IL-17C in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: IL17C Recombinant monoclonal antibody, PBS Only (Capture) 87706-2-PBS (251249G4). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL17C Recombinant monoclonal antibody, PBS Only (Detector) 87706-1-PBS (251249F9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP03144-1, Human IL-17C Recombinant Matched Antibody Pair - PBS only. 87706-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2788. 87706-1-PBS was HRP conjugated as the detection antibody. Range: 15.6-1000 pg/mL



HCT-116 cells were cultured unstimulated or stimulated with 15 ng/mL IL-1 β for 1 day. The mean IL-17C concentration was determined to be 10.2 pg/mL in unstimulated cell culture supernatant, 453.2 pg/mL in stimulated cell culture supernatant.