

For Research Use Only

Mouse CD224/GGT1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03153-2

Capture Antibody Information

Catalog Number:
87659-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
253090F7

Reactivity:
mouse

Immunogen Catalog Number:
Eg6301

Conjugate:
Unconjugated

Full name:
gamma-glutamyltransferase 1

Gene ID:
14598

Detection Antibody Information

Catalog Number:
87659-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
253090A12

Reactivity:
mouse

GenBank:
NM_001305992.1

Immunogen Catalog Number:
Eg6301

Conjugate:
Unconjugated

Full name:
gamma-glutamyltransferase 1

Gene ID:
14598

Applications

Tested Applications:
Sandwich ELISA

Range:
0.39-25 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03153-2 targets CD224/GGT1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Mouse CD224/GGT1 Recombinant monoclonal antibody, PBS Only (Capture) 87659-3-PBS (253090F7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD224/GGT1 Recombinant monoclonal antibody, PBS Only (Detector) 87659-2-PBS (253090A12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

