

For Research Use Only

Ataxin-3/ATXN3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03223-1

Capture Antibody Information

Catalog Number:
87817-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
253467B1

Reactivity:
human

Immunogen Catalog Number:
Ag41832

Conjugate:
Unconjugated

Full name:
ataxin 3

Gene ID:
4287

Detection Antibody Information

Catalog Number:
87817-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
253467E9

Reactivity:
human

GenBank:
Immunogen Catalog Number:
Ag41832

Conjugate:
Unconjugated

Full name:
ataxin 3

Gene ID:
4287

Applications

Tested Applications:
Cytometric bead array

Range:
1.562-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03223-1 targets Ataxin-3/ATXN3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MU154079-4 ATXN3 Recombinant monoclonal antibody, PBS Only (Capture) 87817-3-PBS (253467B1). 100 µg. Concentration 1 mg/ml.

Detection antibody: MU154079-4 ATXN3 Recombinant monoclonal antibody, PBS Only (Detector) 87817-2-PBS (253467E9). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03223-1, Ataxin-3/ATXN3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87817-3-PBS. Detection antibody: 87817-2-PBS. Standard: Ag41832. Range: 1.562-100 ng/mL



The mean Ataxin-3/ATXN3 concentration was determined to be 567.7 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load, 353.4 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load, 132.2 ng/mL in U-251 cell extract based on a 1.2 mg/mL extract load.