

For Research Use Only

# Gephyrin/GPHN Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03251-1

## Capture Antibody Information

Catalog Number:  
87903-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
260132C5

Reactivity:  
human

Immunogen Catalog Number:  
Ag3370

Conjugate:  
Unconjugated

Full name:  
gephyrin

Gene ID:  
10243

## Detection Antibody Information

Catalog Number:  
87903-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
260132A2

Reactivity:  
human

GenBank:  
BC030016

Immunogen Catalog Number:  
Ag3370

Conjugate:  
Unconjugated

Full name:  
gephyrin

Gene ID:  
10243

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP03251-1 targets Gephyrin/GPHN in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GPHN Recombinant monoclonal antibody, PBS Only (Capture) 87903-4-PBS (260132C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: GPHN Recombinant monoclonal antibody, PBS Only (Detector) 87903-3-PBS (260132A2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP03251-1, Gephyrin/GPHN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87903-4-PBS. Detection antibody: 87903-3-PBS. Standard: Ag3370. Range: 0.625-80 ng/mL.



The mean Gephyrin/GPHN concentration was determined to be 79.1 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 58.2 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 23.1 ng/mL in Caco-2 cell extract based on a 1.5 mg/mL extract load.