

For Research Use Only



NFATC2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP03256-1

Capture Antibody Information

Catalog Number:
87895-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
253473B2
Reactivity:
human
Immunogen Catalog Number:
Ag16849

Conjugate:
Unconjugated
Full name:
nuclear factor of activated T-cells,
cytoplasmic, calcineurin-dependent 2
Gene ID:
4773

Detection Antibody Information

Catalog Number:
87895-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
253473H4
Reactivity:
human
GenBank:
BC136418
Immunogen Catalog Number:
Ag16849

Conjugate:
Unconjugated
Full name:
nuclear factor of activated T-cells,
cytoplasmic, calcineurin-dependent 2
Gene ID:
4773

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03256-1 targets NFATC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFATC2 Recombinant monoclonal antibody, PBS Only (Capture) 87895-3-PBS (253473B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: NFATC2 Recombinant monoclonal antibody, PBS Only (Detector) 87895-2-PBS (253473H4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03256-1, NFATC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87895-3-PBS. Detection antibody: 87895-2-PBS. Standard: Ag16849. Range: 0.156-20 ng/mL



The mean NFATC2 concentration was determined to be 2.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 27.7 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, RAW267.4 cell extract were served as a negative control.