

E2F1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03288-1**

Capture Antibody Information

Catalog Number:

87313-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252660F3

Reactivity:

human

Immunogen Catalog Number:

Ag17363

Conjugate:

Unconjugated

Full name:

E2F transcription factor 1

Gene ID:

1869

Detection Antibody Information

Catalog Number:

87313-7-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

252660G9

Reactivity:

human

GenBank:

BC050369

Immunogen Catalog Number:

Ag17363

Conjugate:

Unconjugated

Full name:

E2F transcription factor 1

Gene ID:

1869

Applications

Tested Applications:

Sandwich ELISA

Range:

31.3-1000 pg/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03288-1 targets E2F1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: E2F1 Recombinant monoclonal antibody, PBS Only (Capture) 87313-1-PBS (252660F3). 100 µg. Concentration 1 mg/mL.

Detection antibody: E2F1 Recombinant monoclonal antibody, PBS Only (Detector) 87313-7-PBS (252660G9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP03288-1, Human E2F1 Recombinant Matched Antibody Pair - PBS only. 87313-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag17363. 87313-7-PBS was HRP conjugated as the detection antibody. Range: 31.3-1000 pg/mL.



The mean E2F1 concentration was determined to be 1.7 ng/mL in HEK-293 cell extract based on a 1.5 mg/mL extract load, 0.9 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load, 0.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 1.7 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load and 1.2 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.