

For Research Use Only



FASN Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP03330-2

Capture Antibody Information

Catalog Number:
81079-5-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
260308B6

Reactivity:
human

Immunogen Catalog Number:
Ag0975

Conjugate:
Unconjugated

Full name:
fatty acid synthase

Gene ID:
2194

Detection Antibody Information

Catalog Number:
81079-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
260308C7

Reactivity:
human

GenBank:
BC007909

Immunogen Catalog Number:
Ag0975

Conjugate:
Unconjugated

Full name:
fatty acid synthase

Gene ID:
2194

Applications

Tested Applications:
Sandwich ELISA

Range:
0.391-25 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03330-2 targets FASN in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: FASN Recombinant monoclonal antibody, PBS Only (Capture) 81079-5-PBS (260308B6). 100 µg. Concentration 1 mg/mL.

Detection antibody: FASN Recombinant monoclonal antibody, PBS Only (Detector) 81079-3-PBS (260308C7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP03330-2, Human FASN Recombinant Matched Antibody Pair - PBS only. 81079-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag0975. 81079-3-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL

The mean FASN concentration was determined to be 15.53 ng/mL in HEK-293T cell extract based on a 1.20 mg/mL extract load.