

For Research Use Only

# MZB1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP03414-1

## Capture Antibody Information

Catalog Number:  
88201-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
260678B2

Reactivity:  
human

Immunogen Catalog Number:  
Ag42451

Conjugate:  
Unconjugated

Full name:  
hypothetical protein MGC29506

Gene ID:  
51237

## Detection Antibody Information

Catalog Number:  
88201-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
260678E6

Reactivity:  
human

GenBank:  
BC021275

Immunogen Catalog Number:  
Ag42451

Conjugate:  
Unconjugated

Full name:  
hypothetical protein MGC29506

Gene ID:  
51237

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP03414-1 targets MZB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MZB1 Recombinant monoclonal antibody, PBS Only (Capture) 88201-3-PBS (260678B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: MZB1 Recombinant monoclonal antibody, PBS Only (Detector) 88201-2-PBS (260678E6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only



For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP03414-1, MZB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88201-3-PBS. Detection antibody: 88201-2-PBS. Standard: Ag42451. Range: 0.625-80 ng/mL

The mean MZB1 concentration was determined to be 11.6 ng/mL in Raji cell extract based on a 1.0 mg/mL extract load, 25.4 ng/mL in Ramos cell extract based on a 1.0 mg/mL extract load, 2.8 ng/mL in HeLa cell extract based on a 1.0 mg/mL extract load.