

CYP4F11 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03443-1**

Capture Antibody Information

Catalog Number:
88260-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
260726D9
Reactivity:
human
Immunogen Catalog Number:
Ag2906

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 4, subfamily F, polypeptide 11
Gene ID:
57834

Detection Antibody Information

Catalog Number:
88260-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
260726B10
Reactivity:
human
GenBank:
BC016853
Immunogen Catalog Number:
Ag2906

Conjugate:
Unconjugated
Full name:
cytochrome P450, family 4, subfamily F, polypeptide 11
Gene ID:
57834

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03443-1 targets CYP4F11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CYP4F11 Recombinant monoclonal antibody, PBS Only (Capture) 88260-3-PBS (260726D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: CYP4F11 Recombinant monoclonal antibody, PBS Only (Detector) 88260-2-PBS (260726B10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP03443-1, CYP4F11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88260-3-PBS. Detection antibody: 88260-2-PBS. Standard: Ag2906. Range: 0.312-40 ng/mL



The mean CYP4F11 concentration was determined to be 3.0 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load, Ramos cell extract were served as a negative control.