

Syntabulin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP03456-1

Capture Antibody Information

Catalog Number:
88360-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
260730A9

Reactivity:
human

Immunogen Catalog Number:
Ag9344

Conjugate:
Unconjugated

Full name:
Golgi-localized protein

Gene ID:
55638

Detection Antibody Information

Catalog Number:
88360-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
260730D4

Reactivity:
human

GenBank:
BC012082

Immunogen Catalog Number:
Ag9344

Conjugate:
Unconjugated

Full name:
Golgi-localized protein

Gene ID:
55638

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03456-1 targets Syntabulin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GOLSYN Recombinant monoclonal antibody, PBS Only (Capture) 88360-2-PBS (260730A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: GOLSYN Recombinant monoclonal antibody, PBS Only (Detector) 88360-1-PBS (260730D4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03456-1, Syntabulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88360-2-PBS. Detection antibody: 88360-1-PBS. Standard: Ag9344. Range: 0.312-40 ng/mL



The mean Syntabulin concentration was determined to be 4.8 ng/mL in HT-29 cell extract based on a 1.2 mg/mL extract load, 2.5 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, HL-60 cell extract were served as a negative control.