

# MAGOH Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03492-1**

## Capture Antibody Information

**Catalog Number:**  
88341-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
260941F11  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag3004

**Conjugate:**  
Unconjugated  
**Full name:**  
mago-nashi homolog, proliferation-associated (Drosophila)  
**Gene ID:**  
4116

## Detection Antibody Information

**Catalog Number:**  
88341-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
260941B1  
**Reactivity:**  
human  
**GenBank:**  
BC018211  
**Immunogen Catalog Number:**  
Ag3004

**Conjugate:**  
Unconjugated  
**Full name:**  
mago-nashi homolog, proliferation-associated (Drosophila)  
**Gene ID:**  
4116

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.938-60 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP03492-1 targets MAGOH in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MAGOH Recombinant monoclonal antibody, PBS Only (Capture) 88341-1-PBS (260941F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: MAGOH Recombinant monoclonal antibody, PBS Only (Detector) 88341-2-PBS (260941B1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com) W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP03492-1, MAGOH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88341-1-PBS. Detection antibody: 88341-2-PBS. Standard: Ag3004. Range: 0.938-60 ng/mL



The mean MAGOH concentration was determined to be 93.0 ng/mL in HL-60 cell extract based on a 1.0 mg/mL extract load, 117.5 ng/mL in THP-1 cell extract based on a 1.0 mg/mL extract load, 35.0 ng/mL in A549 cell extract based on a 1.0 mg/mL extract load.