

Centrin 2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03517-2**

Capture Antibody Information

Catalog Number:

88406-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

260861E4

Reactivity:

human

Immunogen Catalog Number:

Ag8638

Conjugate:

Unconjugated

Full name:

centrin, EF-hand protein, 2

Gene ID:

1069

Detection Antibody Information

Catalog Number:

88406-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

260861G6

Reactivity:

human

GenBank:

BC005334

Immunogen Catalog Number:

Ag8638

Conjugate:

Unconjugated

Full name:

centrin, EF-hand protein, 2

Gene ID:

1069

Applications

Tested Applications:

Sandwich ELISA

Range:

39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP03517-2 targets Centrin 2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Centrin 2 Recombinant monoclonal antibody, PBS Only (Capture) 88406-4-PBS (260861E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: Centrin 2 Recombinant monoclonal antibody, PBS Only (Detector) 88406-2-PBS (260861G6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP03517-2, Human Centrin 2 Recombinant Matched Antibody Pair - PBS only. 88406-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag8638. 88406-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



The mean Centrin 2 concentration was determined to be 8.07 ng/mL in U-87 MG cell extract based on a 1.2 mg/mL extract load and 1.29 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load.