

BLVRA Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP03540-1**

Capture Antibody Information

Catalog Number: 88427-3-PBS	Clone ID: 260945C3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: biliverdin reductase A
Isotype: IgG	Immunogen Catalog Number: Ag1166	Gene ID: 644
Purification Method: Protein A purification		

Detection Antibody Information

Catalog Number: 88427-2-PBS	Clone ID: 260945A12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: biliverdin reductase A
Isotype: IgG	GenBank: BC008456	Gene ID: 644
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1166	

Applications

Tested Applications: Cytometric bead array	Range: 0.469-60 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP03540-1 targets BLVRA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BLVRA Recombinant monoclonal antibody, PBS Only (Capture) 88427-3-PBS (260945C3). 100 µg. Concentration 1 mg/mL.

Detection antibody: BLVRA Recombinant monoclonal antibody, PBS Only (Detector) 88427-2-PBS (260945A12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP03540-1, BLVRA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88427-3-PBS. Detection antibody: 88427-2-PBS. Standard: Ag1166. Range: 0.469-60 ng/mL.



The mean BLVRA concentration was determined to be 48.7 ng/mL in U2OS cell extract based on a 1.2 mg/mL extract load, 100.6 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, Ramos cell extract were served as a negative control.