

For Research Use Only

IL-17F Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50006-1

Capture Antibody Information

Catalog Number:
68666-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2C7D11

Reactivity:
human

Immunogen Catalog Number:
HZ-1116

Conjugate:
Unconjugated

Full name:
interleukin 17F

Gene ID:
112744

Detection Antibody Information

Catalog Number:
68666-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1C5C6

Reactivity:
human

GenBank:
Immunogen Catalog Number:

HZ-1116

Conjugate:
Unconjugated

Full name:
interleukin 17F

Gene ID:
112744

Applications

Tested Applications:
Sandwich ELISA

Range:
62.5-4000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50006-1 targets IL-17F in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: IL-17F Monoclonal antibody, PBS Only (Capture) 68666-1-PBS (2C7D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-17F Monoclonal antibody, PBS Only (Detector) 68666-2-PBS (1C5C6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50006-1, IL-17F Monoclonal Matched Antibody Pair, PBS Only. 68666-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard HZ-1116. 68666-2-PBS was HRP conjugated as the detection antibody Range: 62.5-4000 pg/mL