

KIAA0284 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50021-1**

Capture Antibody Information

Catalog Number:
66749-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G4E3
Reactivity:
Human
Immunogen Catalog Number:
Ag26852

Conjugate:
Unconjugated
Full name:
KIAA0284
Gene ID:
283638

Detection Antibody Information

Catalog Number:
66749-1-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
1D4B8
Reactivity:
Human
GenBank:
BC112928
Immunogen Catalog Number:
Ag26852

Conjugate:
Unconjugated
Full name:
KIAA0284
Gene ID:
283638

Applications

Tested Applications:
Sandwich ELISA

Range:
0.313-20 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50021-1 targets KIAA0284 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: KIAA0284 Monoclonal antibody, PBS Only (Capture) 66749-2-PBS (1G4E3). 100 µg. Concentration 1 mg/mL.

Detection antibody: KIAA0284 Monoclonal antibody, PBS Only (Detector) 66749-1-PBS (1D4B8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve ofMP50021-1, KIAA0284 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66749-2-PBS. Detection antibody: HRP-conjugate66749-1-PBS. Standard: Ag26852. Range: 0.313-20 ng/mL.