

PD-L1/CD274 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50090-1**

Capture Antibody Information

Catalog Number:
68745-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C12F12
Reactivity:
human

Conjugate:
Unconjugated
Full name:
CD274 molecule
Gene ID:
29126

Detection Antibody Information

Catalog Number:
68745-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2D4A12
Reactivity:
human
GenBank:
BC074984

Conjugate:
Unconjugated
Full name:
CD274 molecule
Gene ID:
29126

Applications

Tested Applications:
Sandwich ELISA

Range:
0.049-1.563 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50090-1 targets PD-L1/CD274 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: PD-L1/CD274 Monoclonal antibody, PBS Only (Capture) 68745-1-PBS (2C12F12). 100 µg. Concentration 1 mg/mL.

Detection antibody: PD-L1/CD274 Monoclonal antibody, PBS Only (Detector) 68745-2-PBS (2D4A12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50090-1, PD-L1/CD274 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68745-1-PBS. Detection antibody: HRP-conjugated 68745-2-PBS. Standard:Eg31414. Range:0.049-1.563 ng/mL