

BCAS2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50169-1**

Capture Antibody Information

Catalog Number:
68003-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2G5B11
Reactivity:
Human
Immunogen Catalog Number:
Ag21309

Conjugate:
Unconjugated
Full name:
breast carcinoma amplified sequence 2
Gene ID:
10286

Detection Antibody Information

Catalog Number:
68003-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1C1E2
Reactivity:
Human
GenBank:
BC005285
Immunogen Catalog Number:
Ag21309

Conjugate:
Unconjugated
Full name:
breast carcinoma amplified sequence 2
Gene ID:
10286

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50169-1 targets BCAS2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCAS2 Monoclonal antibody, PBS Only (Capture) 68003-2-PBS (2G5B11). 100 µg. Concentration 1 mg/mL.

Detection antibody: BCAS2 Monoclonal antibody, PBS Only (Detector) 68003-3-PBS (1C1E2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50169-1, BCAS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68003-2-PBS. Detection antibody: 68003-3-PBS. Standard:Ag21309. Range: 3.125-100 ng/mL.