

For Research Use Only

AFF1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50192-1

Capture Antibody Information

Catalog Number:
68817-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C2G3
Reactivity:
Human
Immunogen Catalog Number:
Ag34482

Conjugate:
Unconjugated
Full name:
AF4/FMR2 family, member 1
Gene ID:
4299

Detection Antibody Information

Catalog Number:
68817-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2D6F3
Reactivity:
Human
GenBank:
NM_001166693
Immunogen Catalog Number:
Ag34482

Conjugate:
Unconjugated
Full name:
AF4/FMR2 family, member 1
Gene ID:
4299

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50192-1 targets AFF1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AFF1 Monoclonal antibody, PBS Only (Capture) 68817-1-PBS (1C2G3). 100 µg. Concentration 1 mg/mL.

Detection antibody: AFF1 Monoclonal antibody, PBS Only (Detector) 68817-2-PBS (2D6F3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50192-1, AFF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68817-1-PBS. Detection antibody: 68817-2-PBS. Standard:Ag34482. Range: 0.781-100 ng/mL.