

# AMACR/p504S Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50263-1**

## Capture Antibody Information

Catalog Number:  
60240-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein A purification

Clone ID:  
5F4D7

Reactivity:  
human

Immunogen Catalog Number:  
Ag8720

Conjugate:  
Unconjugated

Full name:  
alpha-methylacyl-CoA racemase

Gene ID:  
23600

## Detection Antibody Information

Catalog Number:  
60240-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein A purification

Clone ID:  
5F5D6

Reactivity:  
human

GenBank:  
BC009471

Immunogen Catalog Number:  
Ag8720

Conjugate:  
Unconjugated

Full name:  
alpha-methylacyl-CoA racemase

Gene ID:  
23600

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50263-1 targets AMACR/p504S in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AMACR/p504S Monoclonal antibody, PBS Only (Capture) 60240-2-PBS (5F4D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: AMACR/p504S Monoclonal antibody, PBS Only (Capture/Detector) 60240-3-PBS (5F5D6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP50263-1, p504S,AMACR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60240-2-PBS. Detection antibody: 60240-3-PBS. Standard:Ag8720. Range: 1.563-200 ng/mL