

TSLP Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50325-2**

Capture Antibody Information

Catalog Number:
68911-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H1G7

Reactivity:
human

Immunogen Catalog Number:
Ag4771

Conjugate:
Unconjugated

Full name:
thymic stromal lymphopoietin

Gene ID:
85480

Detection Antibody Information

Catalog Number:
68911-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E8B1

Reactivity:
human

GenBank:
BC016720

Immunogen Catalog Number:
Ag4771

Conjugate:
Unconjugated

Full name:
thymic stromal lymphopoietin

Gene ID:
85480

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50325-2 targets TSLP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSLP Monoclonal antibody, PBS Only (Capture) 68911-1-PBS (1H1G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: TSLP Monoclonal antibody, PBS Only (Detector) 68911-3-PBS (1E8B1). 100 µg. Concentration 1 mg/mL.

Alternative TSLP matched antibody pairs: MP01287-1, MP01287-2, MP01287-3, MP01287-4, MP50325-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50325-2, TSLP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68911-1-PBS. Detection antibody: 68911-3-PBS. Standard:Ag4771. Range: 0.098-50 ng/mL