

Epiregulin Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50352-1**

Capture Antibody Information

Catalog Number:
68929-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2F6G4

Reactivity:
human

Immunogen Catalog Number:
Ag18915

Conjugate:
Unconjugated

Full name:
epiregulin

Gene ID:
2069

Detection Antibody Information

Catalog Number:
68929-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2B4G4

Reactivity:
human

GenBank:
BC136404

Immunogen Catalog Number:
Ag18915

Conjugate:
Unconjugated

Full name:
epiregulin

Gene ID:
2069

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50352-1 targets Epiregulin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Epiregulin Monoclonal antibody, PBS Only (Capture) 68929-1-PBS (2F6G4). 100 µg. Concentration 1 mg/mL.

Detection antibody: Epiregulin Monoclonal antibody, PBS Only (Capture/Detector) 68929-2-PBS (2B4G4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50352-1, EREG/Epiregulin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68929-1-PBS. Detection antibody: 68929-2-PBS. Standard:Ag18915. Range: 1.563-200 ng/mL.