

CD62L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50444-2**

Capture Antibody Information

Catalog Number:

67444-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

3C3B8

Reactivity:

human

Immunogen Catalog Number:

Ag29297

Conjugate:

Unconjugated

Full name:

selectin L

Gene ID:

6402

Detection Antibody Information

Catalog Number:

67444-4-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1A4B9

Reactivity:

human

GenBank:

BC020758

Immunogen Catalog Number:

Ag29297

Conjugate:

Unconjugated

Full name:

selectin L

Gene ID:

6402

Applications

Tested Applications:

Cytometric bead array

Range:

0.098-6.25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50444-2 targets CD62L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD62L Monoclonal antibody, PBS Only (Capture) 67444-2-PBS (3C3B8). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD62L Monoclonal antibody, PBS Only (Detector) 67444-4-PBS (1A4B9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50444-2, CD62L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67444-2-PBS. Detection antibody: 67444-4-PBS. Standard:Ag29297. Range: 0.098-6.25 ng/mL