

PROK1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50534-3**

Capture Antibody Information

Catalog Number:

60417-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1E12G3

Reactivity:

human

Immunogen Catalog Number:

Ag8079

Conjugate:

Unconjugated

Full name:

prokineticin 1

Gene ID:

84432

Detection Antibody Information

Catalog Number:

60417-4-PBS

Host:

Mouse

Isotype:

IgG2a

Purification Method:

Protein A Magarose purification

Clone ID:

2E10H11

Reactivity:

human

GenBank:

BC025399

Immunogen Catalog Number:

Ag8079

Conjugate:

Unconjugated

Full name:

prokineticin 1

Gene ID:

84432

Applications

Tested Applications:

Cytometric bead array

Range:

3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50534-3 targets PROK1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PROK1 Monoclonal antibody, PBS Only (Capture) 60417-1-PBS (1E12G3). 100 µg. Concentration 1 mg/mL.

Detection antibody: PROK1 Monoclonal antibody, PBS Only (Detector) 60417-4-PBS (2E10H11). 100 µg. Concentration 1 mg/mL.

Alternative PROK1 matched antibody pairs: MP50534-1, MP50534-2, MP50534-4

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50534-3, PROK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60417-1-PBS. Detection antibody: 60417-4-PBS. Standard:Ag8079. Range: 3.125-100 ng/mL.