

For Research Use Only



SQRDL Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP50545-3

Capture Antibody Information

Catalog Number:

68459-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1F10B12

Reactivity:

human

Immunogen Catalog Number:

Ag10762

Conjugate:

Unconjugated

Full name:

sulfide quinone reductase-like (yeast)

Gene ID:

58472

Detection Antibody Information

Catalog Number:

68459-5-PBS

Host:

Mouse

Isotype:

IgG3

Purification Method:

Protein A Magarose purification

Clone ID:

1H7E8

Reactivity:

human

GenBank:

BC016836

Immunogen Catalog Number:

Ag10762

Conjugate:

Unconjugated

Full name:

sulfide quinone reductase-like (yeast)

Gene ID:

58472

Applications

Tested Applications:

Cytometric bead array

Range:

1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50545-3 targets SQRDL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SQRDL Monoclonal antibody, PBS Only (Capture) 68459-2-PBS (1F10B12). 100 µg. Concentration 1 mg/mL.

Detection antibody: SQRDL Monoclonal antibody, PBS Only (Detector) 68459-5-PBS (1H7E8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50545-3, SQRDL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68459-2-PBS. Detection antibody: 68459-5-PBS. Standard:Ag10762. Range: 1.563-100 ng/mL.