

CCL24 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50599-1**

Capture Antibody Information

Catalog Number:

60443-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

2A4E10

Reactivity:

human

Conjugate:

Unconjugated

Full name:

chemokine (C-C motif) ligand 24

Gene ID:

6369

Detection Antibody Information

Catalog Number:

60443-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

2H8D1

Reactivity:

human

GenBank:

NM_002991

Conjugate:

Unconjugated

Full name:

chemokine (C-C motif) ligand 24

Gene ID:

6369

Applications

Tested Applications:

Cytometric bead array

Range:

0.195-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50599-1 targets CCL24 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCL24 Monoclonal antibody, PBS Only (Capture/Detector) 60443-1-PBS (2A4E10). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCL24 Monoclonal antibody, PBS Only (Detector) 60443-2-PBS (2H8D1). 100 µg. Concentration 1 mg/mL.

Alternative CCL24 matched antibody pairs: MP50599-2, MP50599-3

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50599-1, CCL24 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60443-1-PBS. Detection antibody: 60443-2-PBS. Standard:Eg0395. Range: 0.195-100 ng/mL.