

For Research Use Only

PLXNC1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP50652-1

Capture Antibody Information

Catalog Number:
60476-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
3D5E5

Reactivity:
human

Immunogen Catalog Number:
Ag30179

Conjugate:
Unconjugated

Full name:
plexin C1

Gene ID:
10154

Detection Antibody Information

Catalog Number:
60476-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3G12E4

Reactivity:
human

GenBank:
BC156867

Immunogen Catalog Number:
Ag30179

Conjugate:
Unconjugated

Full name:
plexin C1

Gene ID:
10154

Applications

Tested Applications:
N/A

Range:

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50652-1 targets PLXNC1 in immunoassays as a matched antibody pair. Validated in N/A.

Capture antibody: PLXNC1 Monoclonal antibody, PBS Only (Capture) 60476-1-PBS (3D5E5). 100 µg. Concentration 1 mg/mL.

Detection antibody: PLXNC1 Monoclonal antibody, PBS Only (Detector) 60476-2-PBS (3G12E4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50652-1, PLXNC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60476-1-PBS. Detection antibody: 60476-2-PBS. Standard: Ag30179. Range: 3.125-100 ng/mL.