

DMWD Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50758-2**

Capture Antibody Information

Catalog Number:
60536-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F2A2
Reactivity:
human
Immunogen Catalog Number:
Ag36845

Conjugate:
Unconjugated
Full name:
dystrophia myotonica, WD repeat containing
Gene ID:
1762

Detection Antibody Information

Catalog Number:
60536-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1H2G12
Reactivity:
human
GenBank:
BC019266
Immunogen Catalog Number:
Ag36845

Conjugate:
Unconjugated
Full name:
dystrophia myotonica, WD repeat containing
Gene ID:
1762

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50758-2 targets DMWD in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DMWD Monoclonal antibody, PBS Only (Capture/Detector) 60536-2-PBS (1F2A2). 100 µg. Concentration 1 mg/mL.

Detection antibody: DMWD Monoclonal antibody, PBS Only (Detector) 60536-3-PBS (1H2G12). 100 µg. Concentration 1 mg/mL.

Alternative DMWD matched antibody pairs: MP50758-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50758-2, DMWD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60536-2-PBS. Detection antibody: 60536-3-PBS. Standard:Ag36845. Range: 0.098-25 ng/mL