

# BFAR Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50881-2**

## Capture Antibody Information

**Catalog Number:**

60618-3-PBS

**Host:**

Mouse

**Isotype:**

IgG1

**Purification Method:**

Protein G Magarose purification

**Clone ID:**

2F7B9

**Reactivity:**

human

**Immunogen Catalog Number:**

Ag7755

**Conjugate:**

Unconjugated

**Full name:**

bifunctional apoptosis regulator

**Gene ID:**

51283

## Detection Antibody Information

**Catalog Number:**

60618-2-PBS

**Host:**

Mouse

**Isotype:**

IgG1

**Purification Method:**

Protein G Magarose purification

**Clone ID:**

1C1A5

**Reactivity:**

human

**GenBank:**

BC003054

**Immunogen Catalog Number:**

Ag7755

**Conjugate:**

Unconjugated

**Full name:**

bifunctional apoptosis regulator

**Gene ID:**

51283

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.098-6.25 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP50881-2 targets BFAR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BFAR Monoclonal antibody, PBS Only (Capture) 60618-3-PBS (2F7B9). 100 µg. Concentration 1 mg/mL.

Detection antibody: BFAR Monoclonal antibody, PBS Only (Detector) 60618-2-PBS (1C1A5). 100 µg. Concentration 1 mg/mL.

Alternative BFAR matched antibody pairs: MP50881-1

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage buffer:**


PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)W: [ptglab.com](http://ptglab.com)**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP50881-2, BFAR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60618-3-PBS. Detection antibody: 60618-2-PBS. Standard:Ag7755. Range: 0.098-6.25 ng/mL