

SEPT12 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50944-4**

Capture Antibody Information

Catalog Number:
60661-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2C1C2
Reactivity:
human
Immunogen Catalog Number:
Ag4644

Conjugate:
Unconjugated
Full name:
septin 12
Gene ID:
124404

Detection Antibody Information

Catalog Number:
60661-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F8B12
Reactivity:
human
GenBank:
BC035619
Immunogen Catalog Number:
Ag4644

Conjugate:
Unconjugated
Full name:
septin 12
Gene ID:
124404

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP50944-4 targets SEPT12 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SEPT12 Monoclonal antibody, PBS Only (Capture) 60661-1-PBS (2C1C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: SEPT12 Monoclonal antibody, PBS Only (Capture/Detector) 60661-5-PBS (1F8B12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50944-4, SEPT12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60661-1-PBS. Detection antibody: 60661-5-PBS. Standard:Ag4644. Range: 0.781-100 ng/mL.