

FOLR2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51016-1**

Capture Antibody Information

Catalog Number:
60704-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2E11B7
Reactivity:
human

Conjugate:
Unconjugated
Full name:
folate receptor 2 (fetal)
Gene ID:
2350

Detection Antibody Information

Catalog Number:
60704-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2F2C1
Reactivity:
human
GenBank:
NM_000803.5

Conjugate:
Unconjugated
Full name:
folate receptor 2 (fetal)
Gene ID:
2350

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51016-1 targets FOLR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOLR2 Monoclonal antibody, PBS Only (Capture) 60704-1-PBS (2E11B7). 100 µg. Concentration 1 mg/mL.

Detection antibody: FOLR2 Monoclonal antibody, PBS Only (Detector) 60704-2-PBS (2F2C1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51016-1, FOLR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60704-1-PBS. Detection antibody: 60704-2-PBS. Standard:null. Range: null.